



... LABORATORIES ...  
**HYAMED**  
• GENEVA •



Rue de Lausanne 20BIS, 1201 Genève, Switzerland

info@hyamax.com

Tel : +41 22 900 09 32

[www.hyamax.com](http://www.hyamax.com)

[www.hyamax.com](http://www.hyamax.com)



... LABORATORIES ...

**HYAMED**

GENEVA





## ABOUT HYAMED

“We offer innovative and enhanced safe products to customers and physicians.”

Located in Geneva, Switzerland, Laboratories Hyamed SA is specialized in designing, manufacturing, and distributing implantable and injection medical devices in the fields of plastic and aesthetic surgery, dermatology, and anti-aging.

**Hyamax®**  
Crosslinked Hyaluronic Acid Filler

The patented technology of Hyamax® injectable crosslinked hyaluronic acid filler has been developed by scientists and engineers specialized in biopolymers and injectable medical devices for over 15 years, who are from Germany, Switzerland, France, the USA, Canada, and Japan.

The quality management system of Laboratories Hyamed SA is certified against ISO 13485 and Hyamax® has obtained CE certification of EU.

Originated from Switzerland, Hyamed's new product -- monophasic crosslinked sodium hyaluronate gel, applies a world-leading technique.

Hyamed strictly abides by the regulations of the CFDA and the EU and strives for more innovative and safer products.

**Hyamax®**  
Crosslinked Hyaluronic Acid Filler

Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Fine

Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Deep

Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Lips

Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep

Hyamax®  
Crosslinked Hyaluronic Acid Filler Deep

Hyamax®  
Crosslinked Hyaluronic Acid Filler Lips

Hyamax®  
Crosslinked Hyaluronic Acid Filler Pure

Hyamax®  
Crosslinked Hyaluronic Acid Filler Contour

Hyamax®  
Crosslinked Hyaluronic Acid Filler Extra Deep

Hyamax®  
Crosslinked Hyaluronic Acid Filler Volumizer





Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Fine



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Deep



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Fine



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Lips



Hyamax®  
Crosslinked Hyaluronic Acid Filler Plus Extra Deep





# Hyamax<sup>®</sup>

Crosslinked Hyaluronic Acid Filler Plus Extra Deep

## CE Certified Dermal Filler

Hyamax<sup>®</sup> Plus Extra Deep is used for mid and deep dermis injection for scars repair, deep wrinkles, or facial tissue augmentation.

- Temporal zone/nose bridge/cheek/chin
- 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- 2x1.0ml
- 27G needle





# Hyamax<sup>®</sup>

Crosslinked Hyaluronic Acid Filler **Plus Lips**

## CE Certified Dermal Filler

Hyamax<sup>®</sup> Plus Lips is used for lip defects or volume enhancement.

- ▶ Lips
- ▶ 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- ▶ 2x1.0ml
- ▶ 27G needle





# Hyamax<sup>®</sup>

Crosslinked Hyaluronic Acid Filler **Plus Fine**


## CE Certified Dermal Filler

Hyamax<sup>®</sup> Plus Fine is used for superficial dermis injection to correct fine lines and wrinkles or facial defects.

- Forehead lines/glabella lines/tear troughs
- 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- 2x1.0ml
- 27G needle





... LABORATORIES ...  
**HYAMED**  
GENEVA  



**Hyamax<sup>®</sup>**  
Crosslinked Hyaluronic Acid Filler **Plus Deep**

### CE Certified Dermal Filler

Hyamax<sup>®</sup> Plus Deep is used for a mid and deep facial tissue augmentation to correct moderate and severe facial wrinkles and folds.

- ▶ Forehead lines/ glabella lines/ tear troughs
- ▶ 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- ▶ 2x1.0ml
- ▶ 27G needle



... LABORATORIES ...  
**HYAMED**  
GENEVA  




## Hyamax<sup>®</sup> Crosslinked Hyaluronic Acid Filler **Contour**

### Be What You Want To Be

► Hyamax<sup>®</sup> Contour and Contour L-plus are high viscosity dermal fillers for treatment of body rhytides and lipoatrophy with/without 0.3% by weight of lidocaine hydrochloride.

With 0.3% lidocaine: Contour L-plus    Without 0.3% lidocaine: Contour

- Buttocks/breasts
- 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- 2x2ml/1x10ml



## Hyamax<sup>®</sup> Non Crosslinked Hyaluronic Acid Filler **Pure**

### Moisturizing, Fresh, Radiant

► Hyamax<sup>®</sup> Pure is used to improve the hydration, elasticity and radiance of the skin by filling fine lines.

- USED FOR: Face/neck/back of hands
- 1x2.0ml
- 12.5mg/ml, non Crosslinked Hyaluronic Acid





# Hyamax<sup>®</sup>

Crosslinked Hyaluronic Acid Filler Volumizer

**Time Fades, but Beauty Endures.**

Using the latest Swiss patented crosslinking technology, Hyamax<sup>®</sup> Volumizer, the mono-phasic filler, is a three-dimensional network structure. Excellent cohesion enables no tendency of the filler to disperse and flow in water and body fluids, achieving higher stability and the best volumizing and contouring effect.

- Cheek/chin
- 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- 2×1.0ml/2×2.0ml
- 27G TSK needle





# Hyamax<sup>®</sup>

Crosslinked Hyaluronic Acid Filler Skin Booster

## Boost Your Skin, Boost Your Life.

By directly replenishing crosslinked hyaluronic acid to the deep dermis, Hyamax<sup>®</sup> Skin Booster can improve skin structure, activate fibroblast, stimulate collagen regeneration, restore skin vitality, and constantly moisturize the skin, which can effectively reduce the fine lines and make the skin moisturized and radiant.

- ▶ 20mg/ml, BDDE Crosslinked Hyaluronic Acid
- ▶ 2x1.0ml or 2x2.0ml
- ▶ 27G needle
- ▶ 3 treatment sessions at 30-day interval followed by treatment every 2 months





# Hyamax® Mesotherapy Series

## Hyamax® HA18

- | Composition: Hyaluronic Acid + 18 Amino Acids
- › Microneedle/dermaroller
- › 3 sessions spaced out by 20 days

### 6D of Perfection

- Moisturizing: moisturize the skin
- Smooth: smooths the skin
- Delicate: delicate and smooth
- Radiance: improve skin radiance
- Uniform: brighten dark complexion
- Anti-wrinkle: fade fine lines and dry lines





## Hyamax® Mesotherapy Series


Hyamax® COGN3

**Composition:** Hyaluronic Acid +2 Amino Acids  
Recombinant Human Collagen Type III

- › Microneedle/dermaroller
- › 3 sessions spaced out by 20 days

### 6D of Perfection

- Moisturizing
- Anti-wrinkle
- Tightening
- Repair
- Elasticity
- Plumpness



Hyamax® COGN9

**Composition:** Recombinant Human Collagen Type III

- › Microneedle/dermaroller
- › 3 sessions spaced out by 20 days

### 6D of Perfection

- Moisturizing
- Anti-wrinkle
- Tightening
- Repair
- Elasticity
- Plumpness





# Hyamax® Mesotherapy Series

## Hyamax® PDRN Series

- | Composition: Hyaluronic Acid+PDRN+3 other anti-aging active ingredients
- Microneedle/dermaroller
- 3 sessions spaced out by 20 days

### 6D of Perfection

- Anti-inflammation
- Anti-acne
- Repairing the skin
- Anti-wrinkle
- Radiance
- Tightening



	Hyamax® HA18	Hyamax® PDRN60	Hyamax® COGN3	Hyamax® COGN9
 Indications	Fine lines, hydration, radiance	Hydration, anti-inflammation, anti-acne, radiance	Tightening, anti-wrinkle, elasticity, rejuvenation, hydration	Elasticity, rejuvenation, anti-wrinkle, tightening
 HA Concentration	4 mg/ml	4 mg/ml	2 mg/ml	(None)
 Composition	Hyaluronic acid + 18 amino acids	Hyaluronic acid + PDRN + 3 other anti-aging active ingredients	Hyaluronic acid + recombinant human collagen type III	Recombinant human collagen type III
 Protocol	3 sessions spaced out by 20 days	3 sessions spaced out by 20 days	3 sessions spaced out by 20 days	3 sessions spaced out by 20 days
 Packing	4 ml x 5 vials	4 ml x 5 vials	4 ml x 5 vials	4 ml x 5 vials
				



# Hyamax<sup>®</sup> Hyavital


**Composition:** 32 mg(H-HA)+32 mg(L-HA)/2 ml

## Hyamax<sup>®</sup> Hyavital

**Technology:** Non crosslinked high-molecular-weight hyaluronic acid (H-HA) and non crosslinked low-molecular weight hyaluronic acid (L-HA) are mixed under certain conditions including high temperature, to form strong hydrogen bonds to get crosslinked while without using BDDE.

High Molecular Weight HA

Low Molecular Weight HA



**Indication:** Hydration and Lifting for Facial Sagging/Neck Lines/ Skin Aging/ Skin Laxity

- 2x1.0ml/2x2.0ml
- 3 treatment sessions at 30-day interval followed by treatment every 2 months  
Supplied in 1ml or 2ml glass syringes.
- Reversing Time by 5 Points (Bio Aesthetic Points)






## deoxycholic acid injection

20mg/2mL(10mg/mL)  
5x2.0ml/box



### Indication:

To improve the appearance and profile of moderate to severe fat below the chin (submental fat) by dissolving the fat and thus result in perfect facial contour.



deoxycholic acid powder (supply raw material )





Non animal source










# Product Outline

# Safe & Effective

Product				
Indication	Moderate to Deep of Wrinkles	Moderate to Deep of Wrinkles	Lip Augmentation	Fine Lines
Content	20mg/ml	20mg/ml	20mg/ml	20mg/ml
Degree of Cross Linking	**	***	**	**
Needle	27G	27G	27G	27G
Injection Depth	Mid Dermis	Mid-Deep Dermis	Mid Dermis	Superficial Dermis
Duration	6 months	6 months	6 months	6 months

Product					
Indication	Body Contouring	Fine Lines	Cheek/Chin	Medium Lines	Medium Lines
Content	20mg/ml	12.5mg/ml	20mg/ml	20mg/ml	32mg/ml
Degree of Cross Linking	***	/	****	*	/
Needle	27G	27G	27G	30G	29G
Injection Depth	Deep Dermis	Superficial Dermis	Periosteum	Superficial Dermis	Superficial Dermis
Duration	6 months	3 months	12-18 months	6 months	6 months

Hyamax® HA products are a series of high viscosity, biodegradable, non animal-based, clear HA gel using Hyamed’s patented crosslinking process, which creates a thicker, longer lasting gel that improves the ability to treat deep and more pronounced lines and wrinkles such as nasolabial folds and oral commissures. With less than 0.025 IU/mg of endotoxins and less than 2 ppm of protein\*, Hyamax® HA products exceed the highest standard of purity for HA.

The Hyamax® HA products are completely free of all animal source, thus reducing the risk of an immunogenic reaction. The HA utilized in the Hyamax® product line is naturally integrated into the tissue keeping the skin healthy and natural. Because of the non-permanent material contained in Hyamax® HA products, the clients can reassess their needs throughout the course of their treatment as desired.

- The safety of the Hyamax® HA products is our utmost concern. The following tests have been conducted on Hyamax® HA products.
- 4-week muscle implantation
  - Acute systemic injection test
  - Biocompatibility test
  - Cytotoxicity
  - Guinea Pig maximization sensitization
  - Identification, including pH, osmolality, residual BDDE, sterility
  - Intra-cutaneous reactivity
  - Mutagenicity
  - Protein
  - Reverse mutation assay
  - Tests for interactions with blood Hemolysis test
  - Tests for local effect after implantation
  - Acute toxicity
  - Sub chronic toxicity
  - 52-week biocompatibility & toxicity study

Chemical/Physical Characterization  
A 500 µm microscopy photo reveals the smaller particle size and more uniform dispersion of the particles of Hyamax® vs. other popular HA based fillers. This results in a more homogenous gel that is smoother to inject and easier to mold.

500 µm Hyamax®  
500 µm Restylane®  
500 µm Captique®