

CLARA
Premium Filler



HARD MEDIUM VOLUMA BODY





WHY CHOOSE CLARA ?

Good filling effect

High viscoelasticity, deformation is hard to shift, strong support, maintain good volume.

Long-lasting

The high concentration of 24 mg/ml and CLARA cross-linked technology to guarantee a long-lasting effect, can last longer than 12 months.

Safety

American FDA approve materials with the terminal sterilization good biocompatibility.

Relaxed injection experience

Use mono-phase, soft and natural, no particle foreign body feeling. 0.3% lidocaine mixture to avoid the pain and discomfort.



Concentration	24 mg/ml
Volume	1ml×1
Needle	2×25G,27G
Lidocaine	0.3%
Shelf Life	24 months



CLARA HARD PREMIUM FILLER

CLARA Hard ranks among the most popular dermal fillers for deep wrinkle correction and facial sculpting. It has several admirable properties that help it deliver high patient satisfaction and outstanding results. In addition, although it is developed to provide temporary cosmetic rejuvenation results, it is long-lasting than most HA-based injectable implants.

CLARA Hard is one of the best in the armamentarium of CLARA's line-up of skin fillers. It treats severe wrinkles and folds and can do a variety of other things, making it versatile and highly effective. Moreover, it is a cosmetic rejuvenation treatment performed with this product that does not cause any discomfort or harmful side effects for the patients.

Strengths of CLARA Hard

Zero gel migration

Natural volume

Stable, highly cross-linked structure

Has Hyaluronic Acid as its main ingredient (naturally found in the human body)

Concentration 24 mg/ml

Volume 1ml×1

Needle 2×27G

Lidocaine 0.3%

Shelf Life 24 months



CLARA MEDIUM PREMIUM FILLER

CLARA Medium filler are one of the most popular lip fillers from South Korea, manufactured by BEAUTY REVOLUTION company. CLARA Medium filler contained 24 mg/ml of hyaluronic acid (HA) concentration compared to other lip fillers, which contain just 20 mg/ml of HA concentration.

CLARA Medium can be injected into several different areas of the face to treat creases and wrinkles that have been affected by the natural aging process. Because it is produced from premium quality raw materials of non-animal origin, this filler is exceptionally safe and has a minimal side effect profile. This high-quality injectable gel is based on purified hyaluronic acid. It gives the filler much-needed stability, allowing its injection into areas such as the lips, where the muscles contract a lot.

Strengths of CLARA Medium

Pain-free removal of nasolabial folds

Non-surgical correction of deep facial wrinkles

Inexpensive augmentation of lips, cheeks, etc

Instant lips pouting and enhancement of lip lines

Concentration 20 mg/ml

Volume 1mlx1

Needle 2x27G

Lidocaine 0.3%

Shelf Life 24 months



CLARA VOLUMA PREMIUM FILLER

CLARA Voluma is the latest injectable monophasic cross-linked sodium hyaluronate gel developed by Beauty Revolution which is based on Swiss product development technology from equipment to process and from raw materials to finished products. Beauty Revolution is in strict compliance with the CE medical device certification standards for injectable sodium hyaluronate gel as to ensure the safety and efficacy of the product while differentiating it from other products on the market.

Monophasic products are more cohesive and stable than biphasic products, effectively pushing over the tissues and providing significantly better volume restoration of the face than biphasic products.

Strengths of CLARA Voluma

Zygomaticomalar

Anteromedial

Submalar

chin

cheek

Concentration	24 mg/ml
Volume	10ml×1
Lidocaine	0.3%
Shelf Life	24 months



CLARA BODY PREMIUM FILLER

CLARA BODY is a monophasic dermal filler. It is more solid, and durable compare to biophysic dermal filler.

We provide great support and lift to good viscoelasticity. Besides, it is easy to be pushed through a needle.

CLARA Body consists of modified hyaluronic acid(HA) produced by bacteria, and 0.3% lidocaine, formulated to a concentration of 24mg/ml in a physiologic buffer.

IMPORTANT

No needle comes with this product, you will need an easi-trans and a separate sterile needle.

Strengths of CLARA Body

modeling and outlining of muscle contours

highlighting of the collarbone

calf modeling

contouring of the front abdominal wall

thicken the skin of the arms

improving skin tone and firmness

Quality Assurance

GMP workshop

Standardized production process

Production:

Inova filling machine

Getinge sterilizer

Pluemat semi-automatic filling and
sealing equipment

DPB-250G plate for automatic blister packaging

5T / h purified water machine

3T / h water for injection (WFI) machine

1T / h pure steam generator

Area:

Construction area: 8,600 m²

Production area: 5,000 m²

Cleaning rooms: 960 m²

Logistics center: 1,980 m²

R&D and Quality center: 100 m²

Work flow



Raw materials



Filling




Sterilization



Labeling



Packaging



Delivery

Certificates



R&D Advantage



1

10 R&D professionals

2

Over 10% of turnover invested in R&D center

3

Cooperate with famous universities, hospitals and laboratories in Europe

4

Charging of national research projects

Application




Effect comparison



Global Distribution

The revolutionary products of Beauty Revolution Medical have been exported to over eight countries, including Turkey, Georgia, Iran, Iraq, Azerbaijan, Russia, Kazakhstan, and Ukraine. The unwavering innovation and professionalism have been recognized throughout the world.



WHAT CAN WE OFFER YOU?

Excellent Quality & Professional Service

24 hours online support

OEM, ODM and Customized options are all available

Quick delivery

Quick and professional after sales service

CLARA

Premium Filler



Made by order and license
Beauty Revolution Niederlassung Österreich
Add: Schegargasse 9, 1190 Wien